

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Mouse Monoclonal Anti-Human Myeloperoxidase (MPO) Antibody

Catalog Number: 70356

Lot Number:

Quantity: _____ mg (= _____ ml)

Antibody Concentration (A_{280nm}): _____ mg/ml

Subclass: IgG₁

Product Form: Protein A column chromatography purified immunoglobulin fraction in 0.015M KPO₄, pH = 7.40, containing 0.85% NaCl and 0.1% NaN₃.

Storage: At 2-8°C, Do Not Freeze. Precipitation may occur upon freezing.

Specificity: Human Myeloperoxidase.

Suggested Use: For enzyme conjugation (detection).

Pairing Information: This antibody pairs with mouse monoclonal anti-human MPO (P/N: 70355) in the human myeloperoxidase sandwich immunoassay.