

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** Rabbit Polyclonal Anti-Human PSA, Affinity Purified

**Catalog Number:** 70423 (3029)

**Lot Number:**

**Quantity:** \_\_\_\_\_ mg ( = \_\_\_\_\_ ml)

**Antibody Concentration ( $A_{280nm}$ ):** \_\_\_\_\_ mg/ml

**Product Form:** PSA antigen-affinity column chromatography purified immunoglobulin fraction in 0.015M  $KPO_4$  buffer, pH = 7.40, containing 0.85% NaCl, 0.1%  $NaN_3$ .

**Specificity:** Human Prostate-specific antigen

**Cross-Reactivity:** No cross-reactivity with human albumin, AFP, CEA, Ferritin, hCG, PAP, CA125, CA15-3, or CA19-9.

**Storage:** At 2-8° C. DO NOT FREEZE. Precipitation may occur upon freezing.

**Suggested Use:** For enzyme conjugation (detection).

**Pairing Information:** This antibody pair with mouse monoclonal anti-human PSA (P/N: 70338) in the human **Total PSA** sandwich immunoassay.

**BIOCHECK, INC.**

323 Vintage Park Drive

Foster City, CA 94404

Tel: (650) 573-1968 Fax: (650) 573-1969

---