

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** Mouse Monoclonal Anti-Human Breast Cancer Antigen (CA15-3)-HRP-Conjugate

**Catalog Number:** 70563-HRP

**Lot Number:**

**Quantity:** \_\_\_\_\_ ml

**Antibody Concentration ( $A_{280nm}$ ):** \_\_\_\_\_ mg/ml

**Product Form:** In 0.05 M  $KPO_4$  buffer, 0.85% NaCl, pH=7.20, containing protein stabilizers and preservatives

**Storage:** At 2-8°C. **Do not freeze.**

**Dilution Factor:** Greater than 1:1,000 in the human CA15-3 ELISA assay

**Pairing Information:** This conjugate pairs with CA15-3 MoAb, P/N: 70650, in the human CA15-3 sandwich immunoassay.