

# BIOCHECK, INC.

425 Eccles Avenue  
South San Francisco, CA 94080  
Tel: (650) 573-1968 E-mail: info@biocheckinc.com

---

## **PRODUCT SPECIFICATIONS**

<b>Product Name:</b>	<b>Mouse Monoclonal Anti-Human Lipoprotein-Associated Phospholipase A2 (Lp-PLA2)</b>
<b>Catalog Number:</b>	70746
<b>Lot Number:</b>	
<b>Quantity:</b>	
<b>Antibody Concentration (<math>A_{280nm}</math>):</b>	
<b>Subclass:</b>	IgG <sub>1</sub>
<b>Product Form:</b>	PBS, pH 7.4 with 0.09% sodium azide
<b>Storage:</b>	At 2-8°C. <u>Do Not Freeze</u> . Precipitation may occur upon freezing.
<b>Suggested Use:</b>	For enzyme conjugation ( <b>detection</b> ) in the human Lp-PLA2 immunoassay.
<b>Pairing Information:</b>	This antibody pairs with Mouse Monoclonal Anti-Human Lp-PLA2 Part No. 70745