

BIOCHECK, INC.

425 Eccles Avenue
South San Francisco, CA 94080
Tel: (650) 573-1968 E-mail: info@biocheckinc.com

PRODUCT SPECIFICATIONS

Product Name:	Mouse Monoclonal Anti-Human PD-L1
Catalog Number:	70753
Quantity:	0.1 mg, 0.5 mg, or 1.0 mg
Subclass:	IgG ₁
Product Form:	Lyophilized (Carrier-Free)
Storage:	Store at 2-8°C.
Suggested Use:	For enzyme conjugation (detection) in the PD-L1 immunoassay.
Pairing Information:	This antibody pairs with Mouse Monoclonal Anti-PD-L1 Part No. 70751 or Mouse Monoclonal Anti-PD-L1 Part No. 70752.