

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Mouse Monoclonal Anti-Human GI Cancer Antigen
(CA 19-9)-HRP-Conjugate

Catalog Number: 80009 (CA 19-9-HRP)

Lot Number:

Quantity: _____ ml

Antibody Concentration (A_{280nm}): _____ mg/ml

Product Form: In 0.05 M KPO_4 buffer, 0.85% NaCl, pH=7.20.
Containing protein stabilizers and preservatives

Storage: At 2-8°C. **Do not freeze.**

**Dilution Factor:
assay** Greater than 1:1,000 in the human CA 19-9 ELISA