

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** Goat Anti-Human IgE HRP Conjugate

**Catalog Number:** 80018

**Lot Number:**

**Quantity:** \_\_\_\_\_ml

**Antibody Concentration ( $A_{280nm}$ ):** \_\_\_\_\_mg/ml

**Product Form:** In 0.05 M  $KPO_4$  buffer, 0.85% NaCl, pH = 7.20, containing proteins, preservatives and 50% Glycerol.

**Storage:** At  $-20^{\circ}C$

**Dilution Factor:** Greater than 1:1,000 in the human IgE ELISA assay.