

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Mouse Monoclonal Anti- β -hCG HRP Conjugate

Catalog Number: 80023 (5008-HRP)

Lot Number:

Quantity: _____ ml

Antibody Concentration (A_{280nm}): _____ mg/ml

Product Form: In 0.05 M KPO_4 buffer, 0.085% NaCl, pH = 7.20, containing proteins and preservatives.

Storage: At 2-8°C. DO NOT FREEZE. The conjugate may become inactivated upon freezing.

Dilution Factor: Greater than 1:1,000 in the hCG ELISA assay.

Pairing Information: This conjugate pairs with Alpha hCG MoAb (P/N: 70327) in the human hCG sandwich immunoassay.