

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name:	Sheep Polyclonal Anti-Human β-2 Microglobulin-HRP Conjugate
Catalog Number:	80027
Lot Number:	
Quantity:	_____ ml
Antibody Concentration (A_{280nm}):	_____ mg/ml
Product Form:	In 0.05 M KPO_4 buffer, 0.85% NaCl, pH=7.20, containing proteins and preservatives.
Storage:	At 2-8° C. <u>DO NOT FREEZE</u> . The conjugate may become inactivated upon freezing.
Dilution Factor:	Greater than 1:1,000 in the human β -2 microglobulin ELISA assay.
Pairing Information:	This conjugate pairs with mouse monoclonal anti-human β -2 microglobulin (P/N: 70191) in the human β-2 microglobulin sandwich immunoassay.