

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name:	Goat Polyclonal Anti-Human CRP Antibody- HRP Conjugate
Catalog Number:	80059
Lot Number:	
Quantity:	_____ml
Antibody Concentration (A_{280nm}):	_____mg/ml
Product Form:	In 0.05 M KPO_4 buffer, 0.85% NaCl, pH=7.20, containing proteins and preservatives.
Storage:	At 2-8° C. <u>DO NOT FREEZE</u> . The conjugate may become inactivated upon freezing.
Dilution Factor:	Greater than 1:4,000 in the human CRP ELISA assay.
Pairing Information:	This conjugate pairs with goat polyclonal anti-human CRP (P/N: 70243) or rabbit polyclonal anti-human CRP (P/N: 70648) in the human CRP sandwich immunoassay.