

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: T3 Standard Stock Solution

Catalog Number: 90003

Lot Number:

Specific Immunological Activity: _____ ng/ml

Quantity: _____ ml

Product Form: Prepared in Dimethyl Sulfoxide (DMSO).

Storage: Store at 2-8°C.
Before each use, thaw the vial at room temperature and mix well.