

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Myeloperoxidase (MPO) Standard Stock Solution

Catalog Number: 90108

Lot Number:

Source: Human Neutrophil

Quantity: _____ ml

Specific Immunological Activity: _____ ng/ml

Product Form: In buffer solution containing 1.0% (v/v) ProClin-300 as preservative.

Storage: -20°C or below.