

---

**HIGH SENSITIVITY HEPCIDIN 25-ENZYME IMMUNOASSAY TEST KIT**  
**Catalog Number: BC-1307**

---

A high sensitive enzyme immunoassay for the quantitative measurement of Hepcidin-25 in serum or plasma.

- **Sample** Serum, plasma
- **Test Format** 96-well ELISA
- **Intended Use** RUO
- **CE Mark** Y
- **Range** 1 - 81 ng/mL
- **Sensitivity** 0.153 ng/mL
- **Sample Volume** 20 ul
- **Incubation Time(s)** 60 / 30 / 20 min
- **Storage Conditions** 2° C - 8° C