

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: General Conjugate Diluent For
Antibody-HRP Conjugate Dilution Use (β -2 Microglobulin)

Catalog Number: MB-10030

Lot Number:

Quantity: _____ ml

Product Form: In 0.05M Tris buffer, pH=7.60, containing proteins,
preservatives and red dye.

Storage: At 2-8°C. DO NOT FREEZE.

Suggested Use: For use in the dilution of β -2 Microglobulin antibody-HRP conjugate in
the β -2 Microglobulin ELISA assay.