

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** General Conjugate Diluent For  
Antibody-HRP Conjugate Dilution Use (FSH)

**Catalog Number:** MB-10030

**Lot Number:**

**Quantity:** \_\_\_\_\_ ml

**Product Form:** In 0.05M Tris buffer, pH=7.60, containing proteins,  
preservatives and red dye.

**Storage:** At 2-8°C. DO NOT FREEZE.

**Suggested Use:** For use in the dilution of FSH antibody-HRP conjugate in  
the human FSH ELISA assay.