

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: General Conjugate Diluent For
Antibody-HRP Conjugate Dilution Use (Total PSA)

Catalog Number: MB-10030

Lot Number:

Quantity: _____ml

Product Form: In 0.05M Tris buffer, pH=7.60, containing proteins,
preservatives, red dye, and mouse IgG.

Storage: At 2-8°C. **DO NOT FREEZE.**

Suggested Use: For use in the dilution of PSA MoAb--HRP conjugate in
the ELISA assay.