

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Conjugate Diluent For Antibody-HRP Conjugate Dilution Use (Ferritin)

Catalog Number: MB-10030

Lot Number:

Quantity: _____ ml

Product Form: In 0.05M Tris buffer, pH=7.60, containing animal serum proteins, preservatives and red dye.

Storage: At 2-8°C. **DO NOT FREEZE.**

Suggested Use: For use in the dilution of Ferritin antibody-HRP conjugate in the human Ferritin ELISA assay.