

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: CA 15-3 Sample Diluent

Catalog Number: MB-50008

Lot Number:

Quantity: _____ml

Product Form: In 0.015 M KPO₄ buffer, 0.85% (w/v) NaCl, 0.1% (w/v) NaN₃, 0.25% (w/v) BSA, 0.25% (w/v) PVP, 0.005% (w/v) yellow dye, pH = 7.40

Storage: At 2-8°C. DO NOT FREEZE.

Suggested Use: For use in the dilution of CA 15-3 patient samples in the CA 15-3 ELISA test.